

88734

SEARCH REQUEST FORM

Scientific and Technical Information Center

Access DB#

Requester's Full Name: MOLLY CUPERLEY Examiner #: 59757 Date: 03/11/03
 Art Unit: 1641 Phone Number 30 8-4239 Serial Number: PCT/US 02/08808 10/09/04
 Mail Box and Bldg/Room Location: 7E12 8D15 Results Format Preferred (circle) PAPER DISK E-MAIL

f more than one search is submitted, please prioritize searches in order of need.

 Please provide a detailed statement of the search topic, and describe as specifically as possible the subject matter to be searched. Include the elected species or structures, keywords, synonyms, acronyms, and registry numbers, and combine with the concept or utility of the invention. Define any terms that may have a special meaning. Give examples or relevant citations, authors, etc, if known. Please attach a copy of the cover sheet, pertinent claims, and abstract.

Title of Invention: Activity Based Probe Analysis
 Inventors (please provide full names): ACTIVX BIOSCIENCES, INC.
Matthew P. Patricelli
 Earliest Priority Filing Date: 02/05/01

***For Sequence Searches Only* Please include all pertinent information (parent, child, divisional, or issued patent numbers) along with the appropriate serial number.**

1) Please search for each of the compound fragments circled in claim 11 in combination with each of the terms fluoresceⁿ, rhodamine (see claim 10) and the dyes of claim 28.
 2) Please search for the concept of analyzing protein mixtures using fluorescence detection (see claim 1). This concept might also be phrased as analyzing a proteome (claim 14) using a probe and/or combinatorial chemical library (claim 26).

Terms: proteome, combinatorial library (claim 26), protein analysis, probe, fluoresceⁿ, electrophoreⁿ (claim 2), rhodamine (claim 31), carbonyl-*l*-tryptophan, 6-carboxy-tetra-methyl-rhodamine (claim 32), xanthene, naphthylamine, coumarin, cyanine, ~~metal~~ metal chelate, BODIPY, lanthanide cryptate; erbium, terbium, rutherfordium chelates (claim 28)

POINT OF CONTACT:
 PAUL SCHULWITZ
 TECHNICAL INFO. SPECIALIST
 CM1 6806 TEL (703) 305-1954


STAFF USE ONLY

| | Type of Search | Vendors and cost where applicable |
|--|------------------------|-----------------------------------|
| Searcher: _____ | NA Sequence (#) _____ | STN <u>802.25</u> |
| Searcher Phone #: _____ | AA Sequence (#) _____ | Dialog _____ |
| Searcher Location: _____ | Structure (#) <u>2</u> | Questel/Orbit _____ |
| Date Searcher Picked Up: <u>3/17</u> | Bibliographic _____ | Dr.Link _____ |
| Date Completed: <u>3/18</u> | Litigation _____ | Lexis/Nexis _____ |
| Searcher Prep & Review Time: <u>20</u> | Fulltext <u>X</u> | Sequence Systems _____ |
| Clerical Prep Time: _____ | Patent Family _____ | WWW/Internet _____ |
| Online Time: <u>63</u> | Other _____ | Other (specify) _____ |

=> d que

L1

STR

N~G4
@20 21

Ak @22

Considered -
04/10/03
MZC

VAR G1=C/N

VAR G2=CH2/S/O/20

VAR G3=H/NH2

VAR G4=H/22

NODE ATTRIBUTES:

CONNECT IS E3 RC AT 11

CONNECT IS E3 RC AT 13

CONNECT IS E3 RC AT 14

CONNECT IS E3 RC AT 15

CONNECT IS E1 RC AT 22

DEFAULT MLEVEL IS ATOM

GGCAT IS LOC AT 22

DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:


RSPEC 5

NUMBER OF NODES IS 22

STEREO ATTRIBUTES: NONE

L2 22283 SEA FILE=REGISTRY SSS FUL L1

L12 STR



VPA 40-33/31 U
 NODE ATTRIBUTES:
 DEFAULT MLEVEL IS ATOM
 DEFAULT ECLEVEL IS LIMITED

GRAPH ATTRIBUTES:
 RING(S) ARE ISOLATED OR EMBEDDED
 NUMBER OF NODES IS 32

STEREO ATTRIBUTES: NONE

L14 11 SEA FILE=REGISTRY SSS FUL L12

L18 15654 SEA FILE=HCAPLUS ABB=ON PLU=ON (XANTHENE OR NAPHTHYLAMINE OR
 COUMARIN OR CYANINE OR METAL CHELATE OR BODIPY OR LANTHANIDE
 CRYPT? OR ERBIUM OR TERBIUM OR RUTHENIUM OR RHUTHENIUM) (S)DYE

L20 8 SEA FILE=HCAPLUS ABB=ON PLU=ON L2 AND (L14 OR RHODAMIN?) AND
 FLUORESC? AND L18

=> d ibib abs hitstr hitind 1-8

L20 ANSWER 1 OF 8 HCAPLUS COPYRIGHT 2003 ACS
 ACCESSION NUMBER: 2002:615942 HCAPLUS
 DOCUMENT NUMBER: 137:165832
 TITLE: Activity based probe analysis
 INVENTOR(S): Patricelli, Matthew P.
 PATENT ASSIGNEE(S): Activx Biosciences, Inc., USA
 SOURCE: PCT Int. Appl., 62 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|---------------|------|----------|-----------------|----------|
| WO 2002063271 | A2 | 20020815 | WO 2002-US3808 | 20020205 |
| WO 2002063271 | C1 | 20021024 | | |

W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN,
CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,
GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR,
LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH,
PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TN, TR, TT, TZ,
UA, UG, US, UZ, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU,
TJ, TM
RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH,
CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR,
BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG

PRIORITY APPLN. INFO.: US 2001-266687P P 20010205

OTHER SOURCE(S): MARPAT 137:165832

AB The invention concerns methods and compns. are described for analyzing complex protein mixts. using **fluorescent** activity-based probes. In particular, probes that specifically react with and bind to the active form of one or more target proteins are employed. **Fluorescent** signals obtained from the labeled active target proteins can be related to the presence or amt. of active members of the desired target protein class. The methods and compns. described herein can be used, for example, to provide diagnostic information concerning pathogenic states, in identifying proteins that may act as therapeutic targets, and in drug discovery.

IT 446850-50-6P 446850-53-9P 446850-55-1DP,
reaction with rhodamine green 446850-58-4P
446850-61-9P 446850-64-2P 446850-67-5P
446850-69-7DP, reaction with rhodamine green
446850-71-1P 446850-73-3P 446850-76-6P
446850-79-9P 446850-81-3P

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)
(activity based probe anal.)


RN 446850-50-6 HCAPLUS

CN Adenosine, 5'-deoxy-5'-[[[(1,1-dimethylethoxy)carbonyl]amino]-, 2'(or 3')-[(2-aminoethyl)carbamate] (9CI) (CA INDEX NAME)

CM 1


CRN 446850-49-3
CMF C15 H22 N6 O5

Absolute stereochemistry.



CM 2



CRN 109-58-0
CMF C3 H8 N2 O2



RN 446850-53-9 HCAPLUS
CN Adenosine, 5'-amino-5'-deoxy-, monoester with 9-[2-carboxy-4(or 5)-[[[2-(carboxyamino)ethyl]amino]carbonyl]phenyl]-3,6-bis(dimethylamino)xanthylum inner salt (9CI) (CA INDEX NAME)

CM 1


CRN 446850-52-8
CMF C28 H28 N4 O6
CCI IDS



CM 2

CRN 14365-44-7
CMF C10 H14 N6 O3

Absolute stereochemistry.



RN 446850-55-1 HCAPLUS
CN Adenosine, 5'-amino-5'-deoxy-, 2'(or 3')-[(2-aminoethyl)carbamate] (9CI)


(CA INDEX NAME)

CM 1

CRN 14365-44-7

CMF C10 H14 N6 O3

Absolute stereochemistry.



CM 2

CRN 109-58-0

CMF C3 H8 N2 O2

H₂N-CH₂-CH₂-NH-CO₂H

RN 446850-58-4 HCAPLUS


CN Adenosine, 5'-deoxy-5'-[[4-(fluorosulfonyl)benzoyl]amino]-, 2' (or 3')-[[2-[[[3(or 4)-[3,6-bis(dimethylamino)xanthylum-9-yl]-4(or 3)-carboxyphenyl]carbonyl]amino]ethyl]carbamate], inner salt (9CI) (CA INDEX NAME)

CM 1

CRN 446850-57-3



CMF C17 H17 F N6 O6 S

Absolute stereochemistry.



CM 2

CRN 446850-52-8
 CMF C28 H28 N4 O6
 CCI IDS




RN 446850-61-9 HCAPLUS
 CN Adenosine, 5'-deoxy-5'-[[4-(ethenylsulfonyl)benzoyl]amino]-, 2' (or 3')-[[2-[[[3(or 4)-[3,6-bis(dimethylamino)xanthylum-9-yl]-4(or 3)-carboxyphenyl]carbonyl]amino]ethyl]carbamate], inner salt (9CI) (CA INDEX NAME)

CM 1



CRN 446850-60-8
 CMF C19 H20 N6 O6 S

Absolute stereochemistry.



CM 2

CRN 446850-52-8
 CMF C28 H28 N4 O6
 CCI IDS




RN 446850-64-2 HCAPLUS
 CN Adenosine, 5'-deoxy-5'-[(1-oxo-2-propenyl)amino]-, 2'(or 3')-[[2-[[[3(or 4)-[3,6-bis(dimethylamino)xanthylum-9-yl]-4(or 3)-carboxyphenyl]carbonyl]amino]ethyl]carbamate], inner salt (9CI) (CA INDEX NAME)

CM 1



CRN 446850-63-1
 CMF C13 H16 N6 O4

Absolute stereochemistry.



CM 2

CRN 446850-52-8
 CMF C28 H28 N4 O6
 CCI IDS




RN 446850-67-5 HCAPLUS
 CN Adenosine, 5'-[(2-chloro-1-oxo-2-propenyl)amino]-5'-deoxy-, 2'(or 3')-[[2-[[[3(or 4)-[3,6-bis(dimethylamino)xanthylum-9-yl]-4(or 3)-carboxyphenyl]carbonyl]amino]ethyl]carbamate], inner salt (9CI) (CA INDEX NAME)

CM 1



CRN 446850-66-4
 CMF C13 H15 Cl N6 O4

Absolute stereochemistry.



CM 2

CRN 446850-52-8
 CMF C28 H28 N4 O6
 CCI IDS



RN 446850-69-7 HCAPLUS


CN Adenosine, 5'-deoxy-5'-[[4-(fluorosulfonyl)benzoyl]amino]-, 2' (or 3')-[(2-aminoethyl)carbamate] (9CI) (CA INDEX NAME)

CM 1

CRN 446850-57-3

CMF C17 H17 F N6 O6 S


Absolute stereochemistry.



CM 2

CRN 109-58-0

CMF C3 H8 N2 O2




RN 446850-71-1 HCAPLUS

CN Adenosine, 5'-O-[(4-methoxyphenyl)diphenylmethyl]-, 2' (or 3')-[(2-aminoethyl)carbamate] (9CI) (CA INDEX NAME)

CM 1


CRN 51600-11-4
CMF C30 H29 N5 O5

Absolute stereochemistry.



CM 2



CRN 109-58-0
CMF C3 H8 N2 O2



RN 446850-73-3 HCAPLUS
CN Adenosine, 2' (or 3')-[[2-[[3 (or 4)-[3,6-bis(dimethylamino)xanthylium-9-yl]-4 (or 3)-carboxybenzoyl]amino]ethyl]carbamate], inner salt (9CI) (CA INDEX NAME)

CM 1


CRN 446850-52-8
CMF C28 H28 N4 O6
CCI IDS



CM 2

CRN 58-61-7
CMF C10 H13 N5 O4

Absolute stereochemistry.




RN 446850-76-6 HCAPLUS
CN Adenosine, 2' (or 3') - [[2 - [[3 (or 4) - [3,6-bis(dimethylamino)xanthylum-9-yl] - 4 (or 3) - carboxybenzoyl]amino]ethyl]carbamate] 5' - (chloroacetate), inner salt (9CI) (CA INDEX NAME)

CM 1



CRN 446850-75-5
CMF C12 H14 Cl N5 O5

Absolute stereochemistry.



CM 2



CRN 446850-52-8
CMF C28 H28 N4 O6
CCI IDS



RN 446850-79-9 HCAPLUS
 CN Adenosine, 2' (or 3')-[[2-[[3 (or 4)-[3,6-bis(dimethylamino)xanthylum-9-yl]-4 (or 3)-carboxybenzoyl]amino]ethyl]carbamate] 5'-[4-(fluorosulfonyl)benzoate], inner salt (9CI) (CA INDEX NAME)

CM 1


CRN 446850-52-8
 CMF C28 H28 N4 O6
 CCI IDS



CM 2

CRN 57454-44-1
 CMF C17 H16 F N5 O7 S


Absolute stereochemistry.



RN 446850-81-3 HCAPLUS
 CN Adenosine, 5'-amino-5'-deoxy-, 2' (or 3')-[[2-[[[3 (or 4)-[3,6-bis(dimethylamino)xanthylum-9-yl]-4 (or 3)-carboxyphenyl]carbonyl]amino]ethyl]carbamate], inner salt, trifluoroacetate (salt) (9CI) (CA INDEX NAME)

CM 1

CRN 76-05-1
 CMF C2 H F3 O2





CM 2

CRN 446850-53-9
 CMF C38 H40 N10 O8
 CCI IDS

CM 3

CRN 446850-52-8
 CMF C28 H28 N4 O6
 CCI IDS




CM 4

CRN 14365-44-7

CMF C10 H14 N6 O3

Absolute stereochemistry.



IC ICM G01N
 CC 9-14 (Biochemical Methods)
 Section cross-reference(s): 1, 14
 ST protein sepn electrophoresis synthesis **fluorescent** probe drug
 screening
 IT Capillary electrophoresis
Cyanine dyes
 Diagnosis
 Diffusion
 Drug screening
 Dyes
 Electrophoresis apparatus
Fluorescent substances
 Fluorometry
 Functional groups
 Gel electrophoresis
 Labels
 Mass spectrometry
 Pathogen
 Separation
 (activity based probe anal.)

IT Dyes
(metal chelate; activity based probe anal.)

IT Dyes
(naphthylamine; activity based probe anal.)

IT 189200-71-3DP, Rhodamine green, reaction with adenosine derivs.
446833-62-1P 446833-64-3P 446850-50-6P 446850-53-9P
446850-55-1DP, reaction with rhodamine green
446850-58-4P 446850-61-9P 446850-64-2P
446850-67-5P 446850-69-7DP, reaction with
rhodamine green 446850-71-1P 446850-73-3P
446850-76-6P 446850-79-9P 446850-81-3P
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(activity based probe anal.)

L20 ANSWER 2 OF 8 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2002:90063 HCAPLUS

DOCUMENT NUMBER: 136:163716

TITLE: Labeled peptides, proteins and antibodies and
processes and intermediates useful for their
preparation

INVENTOR(S): Hahn, Klaus M.; Touthkine, Alexei; Muthyala, Rajeev;
Kraynov, Vadim; Bark, Steven J.; Burton, Dennis R.;
Chamberlain, Chester

PATENT ASSIGNEE(S): The Scripps Research Institute, USA

SOURCE: PCT Int. Appl., 158 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 3

PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|---|------|----------|-----------------|------------|
| WO 2002008245 | A2 | 20020131 | WO 2001-US22194 | 20010713 |
| WO 2002008245 | A3 | 20030130 | | |
| W: | | | | |
| AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, | | | | |
| CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, | | | | |
| GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, | | | | |
| LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, | | | | |
| RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, | | | | |
| UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM | | | | |
| RW: | | | | |
| GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, | | | | |
| DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR, BF, | | | | |
| BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG | | | | |
| WO 2002028890 | A1 | 20020411 | WO 2000-US26821 | 20000929 |
| W: | | | | |
| AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, | | | | |
| CR, CU, CZ, DE, DK, DM, DZ, EE, ES, FI, GB, GD, GE, GH, GM, HR, | | | | |
| HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, | | | | |
| LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, RO, RU, | | | | |
| SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, | | | | |
| YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM | | | | |
| RW: | | | | |
| GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, | | | | |
| DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, | | | | |
| CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG | | | | |
| US 2002055133 | A1 | 20020509 | US 2001-839577 | 20010420 |
| PRIORITY APPLN. INFO.: | | | US 2000-218113P | P 20000713 |

WO 2000-US26821 W 20000929
 US 2001-279302P P 20010328
 US 2001-839577 A 20010420

OTHER SOURCE(S): MARPAT 136:163716

AB The invention provides peptide synthons having protected functional groups for attachment of desired moieties (e.g. functional mols. or probes). Also provided are peptide conjugates prep'd. from such synthons, and synthon and conjugate prepn. methods including procedures for identifying optimum probe attachment sites. Biosensors are provided having functional mols. that can locate and bind to specific biomols. within living cells. Biosensors can detect chem. and physiol. changes in those biomols. as living cells are moving, metabolizing and reacting to its environment. Methods are included for detecting GTP activation of a Rho GTPase protein using polypeptide biosensors. When the biosensor binds GTP-activated Rho GTPase protein, an environmentally sensitive dye emits a signal of a different lifetime, intensity or wavelength than when not bound. New fluorophores whose **fluorescence** responds to environmental changes are also provided that have improved detection and attachment properties, and that can be used in living cells, or in vitro.

IT 56-65-5, ATP, biological studies


RL: BSU (Biological study, unclassified); BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

RN 56-65-5 HCAPLUS

CN Adenosine 5'-(tetrahydrogen triphosphate) (8CI, 9CI) (CA INDEX NAME)

Absolute stereochemistry.



IC ICM C07K001-00

CC 9-14 (Biochemical Methods)

Section cross-reference(s): 7, 15, 34, 41

ST labeled peptide protein antibody prepn; biosensor targeting biomol living cell probe; GTP activation Rho GTPase detection polypeptide biosensor; fluorophore **fluorescence** probe environmental change living cell

IT Imaging

(FLAIR (**fluorescent** activation indicator for Rho proteins); labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT Proteins

RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(cyan **fluorescent** protein, conjugates, polypeptide biosensor

contg.; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Fluorescent dyes**

(**cyanine**, conjugates with peptides; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Proteins**

RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(enhanced green **fluorescent** protein, conjugates, polypeptide biosensor contg.; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Proteins**

RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(enhanced yellow green **fluorescent** protein, conjugates, polypeptide biosensor contg.; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Cyanine dyes**

(**fluorescent**, conjugates with peptides; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Fluorescent substances**

(fluorophores, for detecting changes in responses of living cells to environment; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Rho protein (G protein)**

RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(fusion proteins with **fluorescent** proteins; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Proteins**

RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(green **fluorescent**, conjugates, polypeptide biosensor contg.; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT **Biosensors**

Blood serum

Cell

Cell migration

Endoplasmic reticulum

Fibroblast

Fluorescence

Fluorescence excitation

Fluorescence resonance energy transfer

Fluorescent dyes

Genetic vectors

Human
 Phosphorescence
 Phosphorescent substances
 Signal transduction, biological
 Stress, animal
 Stress, microbial
 Stress, plant
 (labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT Fusion proteins (chimeric proteins)
 RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)
 (of Rho GTPase protein and **fluorescent** proteins; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT Proteins
 RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)
 (red **fluorescent** protein, conjugates, polypeptide biosensor contg.; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT Proteins
 RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)
 (yellow green **fluorescent** protein, conjugates, polypeptide biosensor contg.; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT 9059-32-9DP, GTPase, conjugates with **fluorescent** proteins
 RL: ARG (Analytical reagent use); BPN (Biosynthetic preparation); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)
 (GTP-activated Rho; labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT 65-61-2DP, Acridine Orange, conjugates with peptides 1239-45-8DP, Ethidium Bromide, conjugates with peptides 1325-87-7DP, Cascade Blue, conjugates with peptides 1461-15-0DP, Calcein, conjugates with peptides 2321-07-5DP, **Fluorescein**, conjugates with peptides 2768-89-0DP, **Rhodamine X**, conjugates with peptides 3520-42-1DP, Lissamine **Rhodamine B**, conjugates with peptides 7059-24-7DP, Chromomycin A3, conjugates with peptides 7240-37-1DP, 7-AAD, conjugates with peptides 10199-91-4DP, NBD, conjugates with peptides 18378-89-7DP, Mithramycin, conjugates with peptides 23491-45-4DP, Hoechst 33258, conjugates with peptides 23491-52-3DP, Hoechst 33342, conjugates with peptides 25535-16-4DP, Propidium Iodide, conjugates with peptides 30230-57-0DP, conjugates with peptides 41085-99-8DP, conjugates with peptides 43070-85-5DP, Hydroxycoumarin, conjugates with peptides 47165-04-8DP, DAPI, conjugates with peptides 51908-46-4DP, Dansyl aziridine, conjugates with peptides 70281-37-7DP, Tetramethylrhodamine, conjugates with peptides 76421-73-3DP, Monochlorobimane, conjugates with peptides 76433-29-9DP, LDS 751,

conjugates with peptides 82354-19-6DP, Texas Red, conjugates with peptides 82446-52-4DP, Lucifer Yellow, conjugates with peptides 96314-96-4DP, Indo-1, conjugates with peptides 96314-98-6DP, Fura-2, conjugates with peptides 107091-89-4DP, Thiazole Orange, conjugates with peptides 107347-53-5DP, TRITC, conjugates with peptides 112117-57-4DP, conjugates with peptides 123632-39-3DP, Fluo-3, conjugates with peptides 126208-12-6DP, Carboxy-SNARF-1, conjugates with peptides 143245-02-7DP, conjugates with peptides 143413-84-7DP, TOTO-1, conjugates with peptides 143413-85-8DP, YOYO-1, conjugates with peptides 146368-15-2DP, Cy5, conjugates with peptides 146368-16-3DP, Cy3, conjugates with peptides 149838-22-2DP, FM 1-43, conjugates with peptides 153967-04-5DP, SNARF, conjugates with peptides 157199-59-2DP, TO-PRO-1, conjugates with peptides 157199-63-8DP, TO-PRO-3, conjugates with peptides 165599-63-3DP, BODIPY-FL, conjugates with peptides 166196-17-4DP, TOTO-3, conjugates with peptides 169799-14-8DP, Cy7, conjugates with peptides 194100-76-0DP, SYTOX Green, conjugates with peptides 204934-16-7DP, BODIPY TR, conjugates with peptides 237752-36-2DP, Red 613, conjugates with peptides 247145-11-5DP, Alexa-532, conjugates with peptides 287384-28-5DP, BODIPY TMR, conjugates with peptides 324767-53-5DP, SYTOX Orange, conjugates with peptides 396076-95-2DP, TruRed, conjugates with peptides 396077-00-2DP, SYTOX Blue, conjugates with peptides

RL: ARG (Analytical reagent use); BUU (Biological use, unclassified); PRP (Properties); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)
(labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

IT 56-65-5, ATP, biological studies 86-01-1, GTP 22537-22-0, Magnesium ion, biological studies 142805-58-1, MEK kinase
RL: BSU (Biological study, unclassified); BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)
(labeled peptides, proteins and antibodies and processes and intermediates useful for prepn.)

L20 ANSWER 3 OF 8 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2001:561635 HCAPLUS

DOCUMENT NUMBER: 135:149435

TITLE: Coumarin 6, hypericin, resorufins, and flavins:

suitable chromophores for **fluorescence**

correlation spectroscopy of biological molecules

AUTHOR(S): Benes, Martin; Hudecek, Jiri; Anzenbacher, Pavel; Hof, Martin

CORPORATE SOURCE: J. Heyrovsky Institute of Physical Chemistry, Center for Complex Molecular Systems and Biomolecules, Academy of Science of the Czech Republic, Prague, CZ-18223/8, Czech Rep.

SOURCE: Collection of Czechoslovak Chemical Communications (2001), 66(6), 855-869
CODEN: CCCCCA; ISSN: 0010-0765

PUBLISHER: Institute of Organic Chemistry and Biochemistry, Academy of Sciences of the Czech Republic

DOCUMENT TYPE: Journal

LANGUAGE: English

AB In this work we show that the **dyes** coumarin 6, hypericin, 7-O-ethylresorufin and resorufin are suitable for **fluorescence** correlation spectroscopy (FCS) and demonstrate the use of these **dyes** in physiol. relevant protein studies. Since

coumarins are metabolized by cytochromes P 450, the binding of coumarin 6 to cytochrome P 450 3A4 was investigated by FCS. Coumarin 6 appears to be a very bright non-covalent cytochrome P 450 label. When titrating cytochrome P 450 3A4 with coumarin 6, the diffusion time of the coumarin 6/cytochrome P 450 3A4 complex increases roughly two-fold at protein concns. higher than 1 .mu.mol l-1, indicating the formation of cytochrome aggregates. FCS of the FMN (FMN) and FAD (FAD) shows that both endogenous dyes undergo photobleaching. Moreover, FAD appears to be present to great extent, as a non-fluorescent intramol. complex. Anal. of the FCS data of the flavoprotein NADPH-cytochrome P 450 oxidoreductase (mol. wt. 76 500) yielded two components. While the slow component corresponds to a globular protein with the mol. wt. about 75 000, the fast component appears to be due to free diffusing FMN and FAD mols. The amt. of free FMN and FAD increases with increasing laser power. At high laser power a complete photodissocn. of FMN and FAD occurs.


IT 146-14-5, FAD

RL: PEP (Physical, engineering or chemical process); PROC (Process)
(coumarin 6, hypericin, resorufins, and flavins: suitable chromophores
for **fluorescence** correlation spectroscopy of biol. mols.)

RN 146-14-5 HCAPLUS

CN Riboflavin 5'-(trihydrogen diphosphate), P'.fwdarw.5' ester with adenosine
(9CI) (CA INDEX NAME)

Absolute stereochemistry.



CC 9-5 (Biochemical Methods)

IT Diffusion

Fluorometry

Photochemical bleaching

Simulation and Modeling, physicochemical

(coumarin 6, hypericin, resorufins, and flavins: suitable chromophores
for **fluorescence** correlation spectroscopy of biol. mols.)

IT Proteins, specific or class

RL: PEP (Physical, engineering or chemical process); PROC (Process)
(globular; coumarin 6, hypericin, resorufins, and flavins: suitable
chromophores for **fluorescence** correlation spectroscopy of
biol. mols.)

IT 548-04-9, Hypericin 635-78-9, Resorufin 989-38-8, **rhodamine**

6G 5725-91-7, o-Ethylresorufin 38215-36-0, coumarin 6

RL: ARG (Analytical reagent use); PEP (Physical, engineering or chemical
process); ANST (Analytical study); PROC (Process); USES (Uses)

(coumarin 6, hypericin, resorufins, and flavins: suitable chromophores

for **fluorescence** correlation spectroscopy of biol. mols.)
 IT 146-14-5, FAD 146-17-8, Flavin mononucleotide 9035-51-2,
 cytochrome P.450, processes 9039-06-9 329736-03-0, cytochrome P 450
 3A4
 RL: PEP (Physical, engineering or chemical process); PROC (Process)
 (coumarin 6, hypericin, resorufins, and flavins: suitable chromophores
 for **fluorescence** correlation spectroscopy of biol. mols.)
 REFERENCE COUNT: 47 THERE ARE 47 CITED REFERENCES AVAILABLE FOR THIS
 RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L20 ANSWER 4 OF 8 HCAPLUS COPYRIGHT 2003 ACS
 ACCESSION NUMBER: 2000:421344 HCAPLUS
 DOCUMENT NUMBER: 133:55658
 TITLE: A heterogeneous assay for pyrophosphate detection
 using **fluorescent** nucleotide triphosphate
 probes
 INVENTOR(S): Williams, John G. K.
 PATENT ASSIGNEE(S): Li-Cor, Inc., USA
 SOURCE: PCT Int. Appl., 45 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 2
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|--|------|----------|-----------------|----------|
| WO 2000036151 | A1 | 20000622 | WO 1999-US29584 | 19991213 |
| W: AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM RW: GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG | | | | |
| US 6232075 | B1 | 20010515 | US 1999-460304 | 19991213 |
| US 6255083 | B1 | 20010703 | US 1999-460303 | 19991213 |
| US 2001018184 | A1 | 20010830 | US 2001-816720 | 20010321 |
| US 2002115076 | A1 | 20020822 | US 2001-859104 | 20010514 |

PRIORITY APPLN. INFO.:
 US 1998-112078P P 19981214
 US 1999-115496P P 19990111
 US 1999-460303 A3 19991213
 US 1999-460304 A1 19991213

AB Nucleotide triphosphate (NTP) probes contg. a fluorophore attached to the .gamma.-phosphate and a quencher moiety sufficiently proximal to the fluorophore moiety for use in pyrophosphate detection assays are disclosed. These probes exhibit distinguishable **fluorescence** characteristics when the fluorophore is attached to the nucleotide through the .gamma.-phosphate and when it is unattached to the nucleotide. The present invention also provides kits and integrated systems/methods for practicing the assays described herein. The method is based on incorporation of the NTP into a nucleic acid primer strand using polymerase immobilized on a solid support, thereby releasing the **fluorescent** probe. A change in **fluorescence** characteristics is detected through either **fluorescent** intensity

or lifetime measurement.

IT **56-65-5D**, Adenosine triphosphate, **fluorescent** labeled, biological studies


RL: BPR (Biological process); BSU (Biological study, unclassified); BIOL (Biological study); PROC (Process)

(a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

RN 56-65-5 HCAPLUS

CN Adenosine 5'-(tetrahydrogen triphosphate) (8CI, 9CI) (CA INDEX NAME)

Absolute stereochemistry.



IC ICM C12Q001-68

ICS C12P019-34; C07H019-00; C07H021-00; C07H021-02; C07H021-04

CC 9-15 (Biochemical Methods)

Section cross-reference(s): 3

ST pyrophosphate detection assay **fluorescent** NTP probe primer polymerase

IT **Fluorescent** probes

Nucleic acid amplification (method)

PCR (polymerase chain reaction)

Test kits

(a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Deoxyribonucleoside triphosphates

Nucleoside triphosphates

RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)

(a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Primers (nucleic acid)

RL: ARG (Analytical reagent use); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); USES (Uses)

(a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Crosslinking agents

(alkylene or alkynylamino, for fluorophores; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Resonant energy transfer

(between fluorophore and quencher via; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT **Cyanine dyes**

(fluorophore; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Plate glass
Polyamide fibers, uses
RL: DEV (Device component use); USES (Uses)
(immobilization of polymerase and DNA on; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT DNA
RL: BPR (Biological process); BSU (Biological study, unclassified); BIOL (Biological study); PROC (Process)
(immobilized; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Enzymes, biological studies
RL: BUU (Biological use, unclassified); CAT (Catalyst use); BIOL (Biological study); USES (Uses)
(immobilized; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT **Fluorescence**
(intensity, measurement of; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Electron transfer
(intramol., between fluorophore and quencher via; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT **Fluorescence** quenching
(intramol., ground-state complex; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT **Fluorescence** quenching
(intramol.; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT **Fluorescence**
(lifetime, measurement of; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Pyrimidine nucleotides
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(linker attachment to 5 position; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Purine nucleotides
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(linker attachment to 7 position; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Glass, uses
RL: DEV (Device component use); USES (Uses)
(porous, immobilization of polymerase and DNA on; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT Phosphate group
(.gamma.-, of NTP, fluorophore attached to; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT 277756-37-3
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(NTP probe; a heterogeneous assay for pyrophosphate detection using **fluorescent** nucleotide triphosphate probes)

IT 14000-31-8, Pyrophosphate

- RL: ANT (Analyte); ANST (Analytical study)
 (a heterogeneous assay for pyrophosphate detection using
fluorescent nucleotide triphosphate probes)
- IT 56-65-5D, Adenosine triphosphate, **fluorescent** labeled,
 biological studies 63-39-8D, Uridine triphosphate, **fluorescent**
 labeled 65-47-4D, CTP, **fluorescent** labeled 86-01-1D,
 Guanosine triphosphate, **fluorescent** labeled 365-08-2D, DTPP,
fluorescent labeled 1927-31-7D, DATP, **fluorescent**
 labeled 2056-98-6D, DCTP, **fluorescent** labeled 2564-35-4D,
 DGTP, **fluorescent** labeled
- RL: BPR (Biological process); BSU (Biological study, unclassified); BIOL
 (Biological study); PROC (Process)
 (a heterogeneous assay for pyrophosphate detection using
fluorescent nucleotide triphosphate probes)
- IT 9012-90-2, DNA polymerase 9014-24-8, DNA dependent RNA polymerase
 9068-38-6, Reverse transcriptase
- RL: BUU (Biological use, unclassified); CAT (Catalyst use); BIOL
 (Biological study); USES (Uses)
 (a heterogeneous assay for pyrophosphate detection using
fluorescent nucleotide triphosphate probes)
- IT 88-68-6, Anthranilamide 91-64-5, Coumarin 7440-27-9D, Terbium,
 chelate, uses 17681-50-4, Reactive Red 4 50402-56-7, EDANS
 76823-03-5, 5-Carboxyfluorescein 138026-71-8, BODIPY
- RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
 (fluorophore; a heterogeneous assay for pyrophosphate detection using
fluorescent nucleotide triphosphate probes)
- IT 79-06-1D, Acrylamide, gel 9003-53-6, Polystyrene 9003-53-6D,
 Polystyrene, avidin coated bead 9004-34-6, Cellulose, uses 9004-54-0,
 Dextran, uses
- RL: DEV (Device component use); USES (Uses)
 (immobilization of polymerase and DNA on; a heterogeneous assay for
 pyrophosphate detection using **fluorescent** nucleotide
 triphosphate probes)
- IT 81-88-9 2321-07-5, **Fluorescein**
- RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
 (quencher or fluorophore; a heterogeneous assay for pyrophosphate
 detection using **fluorescent** nucleotide triphosphate probes)
- IT 99-35-4, Trinitrobenzene 569-64-2, Malachite green 3546-21-2, Ethidium
 6268-49-1 25154-54-5, Dinitrobenzene 25338-56-1, Pyrenebutanoic acid
 70281-37-7, Tetramethyl **rhodamine** 82354-19-6, Texas Red
 202466-51-1
- RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
 (quencher; a heterogeneous assay for pyrophosphate detection using
fluorescent nucleotide triphosphate probes)
- REFERENCE COUNT: 2 THERE ARE 2 CITED REFERENCES AVAILABLE FOR THIS
 RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L20 ANSWER 5 OF 8 HCAPLUS. COPYRIGHT 2003 ACS
 ACCESSION NUMBER: 1999:222928 HCAPLUS
 DOCUMENT NUMBER: 130:264438
 TITLE: Sulfonated xanthene derivatives synthesis and
 applications as **fluorescent** stains
 INVENTOR(S): Mao, Fei; Leung, Wai-Yee; Haugland, Richard P.
 PATENT ASSIGNEE(S): Molecular Probes, Inc., USA
 SOURCE: PCT Int. Appl., 63 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent

LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 2
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|--|------|----------|-----------------|------------|
| WO 9915517 | A1 | 19990401 | WO 1998-US19921 | 19980923 |
| W: AU, CA, JP, US | | | | |
| RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE | | | | |
| US 6130101 | A | 20001010 | US 1997-935963 | 19970923 |
| AU 9895046 | A1 | 19990412 | AU 1998-95046 | 19980923 |
| AU 750380 | B2 | 20020718 | | |
| EP 966458 | A1 | 19991229 | EP 1998-948483 | 19980923 |
| R: AT, BE, CH, DE, DK, ES, FR, GB, IT, LI, NL, SE, IE | | | | |
| JP 2001508494 | T2 | 20010626 | JP 1999-519270 | 19980923 |
| WO 2000017650 | A1 | 20000330 | WO 1999-US22193 | 19990923 |
| W: AU, CA, JP | | | | |
| RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE | | | | |
| AU 9964002 | A1 | 20000410 | AU 1999-64002 | 19990923 |
| PRIORITY APPLN. INFO.: | | | | |
| | | | US 1997-935963 | A 19970923 |
| | | | WO 1998-US19921 | W 19980923 |
| | | | US 1998-209045 | A 19981209 |
| | | | WO 1999-US22193 | W 19990923 |

OTHER SOURCE(S): MARPAT 130:264438


AB The present invention describes **xanthene dyes**, including **rhodamines**, rhodols and **fluoresceins** that are substituted one or more times by a sulfonic acid or a salt of a sulfonic acid. The dyes of the invention, including chem. reactive dyes and dye-conjugates are useful as **fluorescent probes**, particularly in biol. samples.

IT 56-65-5, 5'-ATP, analysis
 RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process); USES (Uses) (sulfonated xanthene derivs. synthesis and applications as **fluorescent stains**)

RN 56-65-5 HCAPLUS

CN Adenosine 5'-(tetrahydrogen triphosphate) (8CI, 9CI) (CA INDEX NAME)

Absolute stereochemistry.



IC ICM C07D311-82

ICS C07D491-14; C07D405-12; C07D491-22; C07H003-06; C07H021-00;

C07H019-04; C07K014-415; G01N001-30

CC 9-15 (Biochemical Methods)
Section cross-reference(s): 6, 27

ST sulfonated **xanthene fluorescent dye** probe
conjugate stain

IT Proteins, specific or class
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(A; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Immunoglobulins
Proteins, specific or class
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(G; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Phycoerythrins
RL: RCT (Reactant); RACT (Reactant or reagent)
(R-phycoerythrins, pyridyldisulfide modified; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Proteins, specific or class
RL: BPR (Biological process); BSU (Biological study, unclassified); SPN (Synthetic preparation); BIOL (Biological study); PREP (Preparation); PROC (Process)
(conjugates, sulfonated xanthene conjugate; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Escherichia coli
(derivatized with amine-reactive sulfonated **xanthene dye**; sulfonated **xanthene** derivs. synthesis and applications as **fluorescent** stains)

IT Staining, biological
Stains, biological
(**fluorescent**; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Gene, animal
RL: ANT (Analyte); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(for actin; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Actins
RL: BSU (Biological study, unclassified); BIOL (Biological study)
(gene for; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Drug delivery systems
(injections, microinjection; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Nerve
(neuron, cell tracing; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Receptors
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(pharmaceutical; sulfonated xanthene derivs. synthesis and applications

as **fluorescent** stains)

IT Organelle
(pinosome; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Artery
(pulmonary, cells; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Animal cell
Bacteria (Eubacteria)
Complexing agents
Drugs
Microparticles
Plant cell
Protista
Virus
Yeast
(sulfonated xanthene conjugate; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Actins
Agglutinins and Lectins
Allophycocyanins
Amino acids, biological studies
Antibodies
Avidins
Biliproteins
Disaccharides
Growth factors, animal
Haptens
Lipids, biological studies
Monosaccharides
Nucleic acids
Nucleotides, biological studies
Peptides, biological studies
Polymers, biological studies
Polysaccharides, biological studies
Toxins
RL: BPR (Biological process); BSU (Biological study, unclassified); SPN (Synthetic preparation); BIOL (Biological study); PREP (Preparation); PROC (Process)
(sulfonated xanthene conjugate; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)

IT Chelating agents
Cytolysis
Drugs
Electroporation
 Fluorescent dyes
 Fluorescent probes
 Fluorescent substances
Ions
Liposomes
Microtubule
Nucleic acid hybridization
Phagocytosis
Staining, biological
Stains, biological
Staphylococcus aureus
Test kits

(sulfonated **xanthene** derivs. synthesis and applications as
fluorescent stains)

IT DNA
RL: ANT (Analyte); ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Actins
Tubulins
mRNA
RL: ANT (Analyte); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Agglutinins and Lectins
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Antibodies
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Antigens
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Avidins
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Carbohydrates, analysis
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Enzymes, analysis
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT Hormone receptors
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as

- fluorescent stains)**
- IT Hormones, animal, analysis
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Peptide receptors
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Peptides, analysis
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Probes (nucleic acid)
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Protein receptors
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Proteins, general, analysis
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT RNA
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Receptors
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Toxins
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent stains)**
- IT Agglutinins and Lectins
RL: RCT (Reactant); RACT (Reactant or reagent)

- (sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT **Dyes**
(**xanthene**; sulfonated **xanthene** derivs. synthesis and applications as **fluorescent** stains)
- IT 178623-12-6DP, **Rhodamine** Red X, conjugates
RL: SPN (Synthetic preparation); PREP (Preparation)
(**Rhodamine** Red X; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 9003-53-6DP, Polystyrene, amine deriv.
RL: ARU (Analytical role, unclassified); BUU (Biological use, unclassified); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)
(**fluorescently** labeled microspheres; sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 58-85-5, Biotin 9013-20-1, Streptavidin 17466-45-4, Phalloidin
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 56-65-5, 5'-ATP, analysis
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); BUU (Biological use, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process); USES (Uses)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 222164-96-7DP, conjugate
RL: ARU (Analytical role, unclassified); BPR (Biological process); BSU (Biological study, unclassified); SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological study); PREP (Preparation); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 222164-86-5P 222164-96-7P
RL: ARU (Analytical role, unclassified); RCT (Reactant); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); RACT (Reactant or reagent)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 222159-90-2P
RL: ARU (Analytical role, unclassified); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 9013-20-1DP, Streptavidin, sulfonated xanthene conjugate
RL: BPR (Biological process); BSU (Biological study, unclassified); SPN (Synthetic preparation); BIOL (Biological study); PREP (Preparation); PROC (Process)
(sulfonated xanthene derivs. synthesis and applications as **fluorescent** stains)
- IT 68-11-1, Mercaptoacetic acid, reactions 117-08-8, Tetrachlorophthalic anhydride 463-71-8, Thiophosgene 552-30-7, Trimellitic anhydride 619-66-9, 4-Carboxybenzaldehyde 652-12-0, Tetrafluorophthalic anhydride 870-46-2, tert-Butyl carbazate 1319-82-0, Aminocaproic acid 5466-84-2, 4-Nitrophthalic anhydride 11032-79-4D, .alpha.-Bungarotoxin, conjugate 35167-99-8D, amino deriv. 37293-51-9, Aminodextran 41175-50-2

51644-96-3 58196-33-1 63095-11-4 93801-18-4D, conjugate
105832-38-0 126695-58-7 163222-21-7D, rhodamine deriv.
conjugate 179898-22-7 220906-39-8 222159-69-5 222159-71-9
222159-75-3 222159-87-7 222164-84-3 222164-97-8

RL: RCT (Reactant); RACT (Reactant or reagent)

(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT 222159-70-8P 222159-72-0P 222159-73-1P 222159-74-2P 222159-79-7P
222159-82-2P 222159-84-4P 222159-85-5P 222164-80-9P 222164-81-0P
222164-92-3P 222164-95-6P 222164-98-9P 222164-99-0P 222165-01-7P
222165-02-8P 222165-04-0P

RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT
(Reactant or reagent)

(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

IT 2321-07-5DP, Fluorescein, conjugates 146397-20-8DP, CY-3,
conjugates 183185-51-5DP, Rhodol Green, conjugates 189200-71-3DP,
Rhodamine Green, conjugates 199745-67-0DP, Texas Red-X,
conjugates 222159-76-4P 222159-78-6P 222159-80-0P 222159-81-1P
222159-82-2DP, conjugate 222159-83-3P 222159-86-6P 222159-92-4DP,
conjugate 222159-93-5DP, conjugate 222164-82-1P 222164-83-2P
222164-86-5DP, conjugate 222164-87-6P 222164-88-7P 222164-91-2P
222164-92-3DP, conjugate 222164-93-4P 222164-95-6DP, conjugate
222165-00-6P 222165-01-7DP, conjugate 222165-04-0DP, spiperone
conjugate

RL: SPN (Synthetic preparation); PREP (Preparation)

(sulfonated xanthene derivs. synthesis and applications as
fluorescent stains)

REFERENCE COUNT: 6 THERE ARE 6 CITED REFERENCES AVAILABLE FOR THIS
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L20 ANSWER 6 OF 8 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1998:186454 HCAPLUS

DOCUMENT NUMBER: 128:227061

TITLE: Alternative dye-labeled primers, ribonucleotides,
deoxyribonucleotides, and dideoxyribonucleotides for
automated DNA analysis and homogeneous
amplification/detection assays

INVENTOR(S): Metzker, Michael L.; Gibbs, Richard A.

PATENT ASSIGNEE(S): Baylor College of Medicine, USA

SOURCE: U.S., 11 pp., Cont.-in-part of U.S. 5,614,386.

CODEN: USXXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 4

PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|--|------|----------|-----------------|----------|
| US 5728529 | A | 19980317 | US 1995-553936 | 19951106 |
| US 5614386 | A | 19970325 | US 1995-494216 | 19950623 |
| WO 9700967 | A1 | 19970109 | WO 1996-US10729 | 19960621 |
| W: AU, CA, JP | | | | |
| RW: AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE | | | | |
| CA 2225531 | AA | 19970109 | CA 1996-2225531 | 19960621 |
| AU 9662886 | A1 | 19970122 | AU 1996-62886 | 19960621 |
| AU 699939 | B2 | 19981217 | | |

EP 833936 A1 19980408 EP 1996-921749 19960621

R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT,
IE, FI

PRIORITY APPLN. INFO.:


US 1995-494216 19950623
US 1995-540228 19951006
US 1995-553936 19951106
US 1996-612036 19960307
WO 1996-US10729 19960621

AB Methods for the use of a class dyes for improved DNA sequencing are provided. A new class of **dyes, BODIPY.RTM.** fluorophores, has been described recently. The parent heterocyclic mol. of the BODIPY.RTM. fluorophores is a dipyrrometheneboron difluoride compd. which is modified to create a broad class of spectrally-discriminating fluorophores. The present invention provides methods for the use of **BODIPY.RTM.** fluorophore-labeled DNA for **dye-primer** sequencing in which the **BODIPY.RTM.s** are attached to the 5' end of sequencing by enzymic incorporation of **fluorescently-labeled** ribonucleotides or deoxyribonucleotides, and provides oligonucleotides labeled with substituted 4,4-difluoro-4-bora-3A,4A-diaza-s-indacene (**BODIPY.RTM.** fluorophore) compds. for performing the Taqman assay. **BODIPY.RTM.** fluorophores have improved spectral characteristics compared to conventional **fluorescein** and **rhodamine dyes**. **BODIPY.RTM.** fluorophores have narrower band width, insensitivity to solvent or pH, and improved photostability; thus, **BODIPY.RTM.** fluorophores lead to improved DNA sequencing and/or detection in any method where electrophoresis and detection of DNA is required. Addnl., the spectral properties of the **BODIPY.RTM.** fluorophores are sufficiently similar in wavelength and intensity to be used with conventional equipment known in the art.


IT **186961-32-ODP**, oligonucleotide primers labeled with **186961-33-IDP**, oligonucleotide primers labeled with
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)
(alternative dye-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

RN 186961-32-0 HCAPLUS

CN Borate(4-), [N-[3-[4-amino-7-[5-O-[hydroxy[[hydroxy(phosphonooxy)phosphinyloxy]phosphinyl]-.beta.-D-ribofuranosyl]-7H-pyrrolo[2,3-d]pyrimidin-5-yl]-2-propynyl]-5-[(5-phenyl-2H-pyrrol-2-ylidene-.kappa.N)methyl]-1H-pyrrole-2-propanamidato(5-)-.kappa.N1]difluoro-, tetrahydrogen, (T-4)- (9CI) (CA INDEX NAME)

● 4 H⁺

RN 186961-33-1 HCAPLUS
 CN Borate(4-), [N-[3-[4-amino-7-[5-O-[hydroxy[[hydroxy(phosphonooxy)phosphiny
 1]oxy]phosphinyl]-.beta.-D-ribofuranosyl]-7H-pyrrolo[2,3-d]pyrimidin-5-yl]-
 2-propynyl]-5-[(3,5-diphenyl-2H-pyrrol-2-ylidene-.kappa.N)methyl]-1H-
 pyrrole-2-propanamidato(5-)-.kappa.N1]difluoro-, tetrahydrogen, (T-4)-
 (9CI) (CA INDEX NAME)

● 4 H⁺

IC ICM C12Q001-68
 ICS C12Q001-70; C07H021-04; C12P019-34
 NCL 435006000
 CC 3-1 (Biochemical Genetics)
 Section cross-reference(s): 33
 IT Deoxyribonucleoside triphosphates
 Nucleoside triphosphates

Primers (nucleic acid)

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(**BODIPY** fluorophore-labeled; alternative **dye**-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

IT **Fluorescent substances**

(**BODIPY**; alternative **dye**-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

IT Deoxyribonucleotides

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(deoxyribodinucleotides, triphosphates, **BODIPY** fluorophore-labeled; alternative **dye**-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

IT 187089-10-7DP, oligonucleotide primers labeled with

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(**BODIPY** 530/550; alternative **dye**-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

IT 150173-78-7DP, oligonucleotide primers labeled with

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(**BODIPY** 576/589; alternative **dye**-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

IT 186961-29-5DP, oligonucleotide primers labeled with

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(**BODIPY** 589/616; alternative **dye**-labeled primers, ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for automated DNA anal. and homogeneous amplification/detection assays)

IT 120718-52-7DP, TAMRA, oligonucleotide primers labeled with

138026-71-8DP, **BODIPY**, oligonucleotide primers labeled with

150152-69-5DP, **BODIPY** 581/591, oligonucleotide primers labeled

with 150173-72-1DP, oligonucleotide primers labeled with

150173-89-0DP, **BODIPY** 564/570, oligonucleotide primers labeled

with 165599-63-3DP, **BODIPY** 503/512, oligonucleotide primers

labeled with 174881-57-3DP, **BODIPY** 523/547, oligonucleotide

primers labeled with 186961-30-8DP, oligonucleotide primers labeled with

186961-31-9DP, oligonucleotide primers labeled with **186961-32-0DP**

, oligonucleotide primers labeled with **186961-33-1DP**,

oligonucleotide primers labeled with 186961-34-2DP, oligonucleotide

primers labeled with 186961-35-3DP, oligonucleotide primers labeled with

186961-36-4DP, oligonucleotide primers labeled with 186961-37-5DP,

oligonucleotide primers labeled with 186961-38-6DP, oligonucleotide

primers labeled with 186961-39-7DP, oligonucleotide primers labeled with

186961-40-0DP, oligonucleotide primers labeled with 186961-41-1DP, \

oligonucleotide primers labeled with 186961-42-2DP, oligonucleotide

primers labeled with 186961-43-3DP, oligonucleotide primers labeled with

186961-44-4DP, oligonucleotide primers labeled with 186961-45-5DP,

oligonucleotide primers labeled with 186961-46-6DP, oligonucleotide

primers labeled with 186961-47-7DP, oligonucleotide primers labeled with

186961-48-8DP, oligonucleotide primers labeled with 186961-49-9DP,
oligonucleotide primers labeled with 186961-50-2DP, oligonucleotide
primers labeled with 186961-51-3DP, oligonucleotide primers labeled with
186961-52-4DP, oligonucleotide primers labeled with 186961-53-5DP,
oligonucleotide primers labeled with

RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(alternative **dye**-labeled primers, ribonucleotides,
deoxyribonucleotides, and dideoxyribonucleotides for automated DNA
anal. and homogeneous amplification/detection assays)

L20 ANSWER 7 OF 8 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1997:172504 HCAPLUS

DOCUMENT NUMBER: 126:167460

TITLE: Alternative dye-labeled primers, ribonucleotides,
deoxyribonucleotides, and dideoxyribonucleotides for
automated DNA analysis and homogeneous
amplification/detection assays

INVENTOR(S): Metzker, Michael L.; Gibbs, Richard A.

PATENT ASSIGNEE(S): Baylor College of Medicine, USA

SOURCE: PCT Int. Appl., 83 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 4

PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|--|------|----------|-----------------|----------|
| WO 9700967 | A1 | 19970109 | WO 1996-US10729 | 19960621 |
| W: AU, CA, JP | | | | |
| RW: AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE | | | | |
| US 5614386 | A | 19970325 | US 1995-494216 | 19950623 |
| US 5861287 | A | 19990119 | US 1995-540228 | 19951006 |
| US 5728529 | A | 19980317 | US 1995-553936 | 19951106 |
| US 5994063 | A | 19991130 | US 1996-612036 | 19960307 |
| AU 9662886 | A1 | 19970122 | AU 1996-62886 | 19960621 |
| AU 699939 | B2 | 19981217 | | |
| EP 833936 | A1 | 19980408 | EP 1996-921749 | 19960621 |
| R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI | | | | |

PRIORITY APPLN. INFO.: US 1995-494216 19950623
US 1995-540228 19951006
US 1995-553936 19951106
US 1996-612036 19960307
WO 1996-US10729 19960621


AB Methods for the use of a class dyes for improved DNA sequencing are
provided. A new class of **dyes, BODIPY.RTM.**
fluorophores, has been described recently. The parent heterocyclic mol.
of the BODIPY.RTM. fluorophores is a dipyrrometheneboron difluoride compd.
which is modified to create a broad class of spectrally-discriminating
fluorophores. The present invention provides methods for the use of
BODIPY.RTM. fluorophore-labeled DNA for **dye**-primer
sequencing in which the **BODIPY.RTM.s** are attached to the 5' end
of sequencing by enzymic incorporation of **fluorescently-labeled**
ribonucleotides or **deoxyribonucleotides**, and provides oligonucleotides
labeled with substituted 4,4-difluoro-4-bora-3A,4A-diaza-s-indacene (

BODIPY.RTM. fluorophore) compds. for performing the Tagman assay. BODIPY.RTM. fluorophores have improved spectral characteristics compared to conventional **fluorescein** and **rhodamine dyes**. BODIPY.RTM. fluorophores have narrower band width, insensitivity to solvent or pH, and improved photostability; thus, BODIPY.RTM. fluorophores lead to improved DNA sequencing and/or detection in any method where electrophoresis and detection of DNA is required. Addnl., the spectral properties of the BODIPY.RTM. fluorophores are sufficiently similar in wavelength and intensity to be used with conventional equipment known in the art.

IT 186961-32-ODP, oligonucleotide primers labeled with
186961-33-IDP, oligonucleotide primers labeled with
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(alternative dye-labeled primers, ribonucleotides,
deoxyribonucleotides, and dideoxyribonucleotides for automated DNA
anal. and homogeneous amplification/detection assays)

RN 186961-32-0 HCAPLUS


CN Borate(4-), [N-[3-[4-amino-7-[5-O-[hydroxy[[hydroxy(phosphonooxy)phosphiny
l]oxy]phosphinyl]-.beta.-D-ribofuranosyl]-7H-pyrrolo[2,3-d]pyrimidin-5-yl]-
2-propynyl]-5-[(5-phenyl-2H-pyrrol-2-ylidene-.kappa.N)methyl]-1H-pyrrole-2-
propanamidato(5-)-.kappa.N1]difluoro-, tetrahydrogen, (T-4)- (9CI) (CA
INDEX NAME)



● 4 H⁺

RN 186961-33-1 HCAPLUS

CN Borate(4-), [N-[3-[4-amino-7-[5-O-[hydroxy[[hydroxy(phosphonooxy)phosphiny
l]oxy]phosphinyl]-.beta.-D-ribofuranosyl]-7H-pyrrolo[2,3-d]pyrimidin-5-yl]-
2-propynyl]-5-[(3,5-diphenyl-2H-pyrrol-2-ylidene-.kappa.N)methyl]-1H-
pyrrole-2-propanamidato(5-)-.kappa.N1]difluoro-, tetrahydrogen, (T-4)-
(9CI) (CA INDEX NAME)



● 4 H⁺

- IC ICM C12P019-34
ICS C12Q001-68; C12Q001-70; C07H019-04
- CC 3-1 (Biochemical Genetics)
Section cross-reference(s): 33
- IT Deoxyribonucleoside triphosphates
Nucleoside triphosphates
Primers (nucleic acid)
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(BODIPY fluorophore-labeled; alternative dye
-labeled primers, ribonucleotides, deoxyribonucleotides, and
dideoxyribonucleotides for automated DNA anal. and homogeneous
amplification/detection assays)
- IT **Fluorescent substances**
(BODIPY; alternative dye-labeled primers,
ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for
automated DNA anal. and homogeneous amplification/detection assays)
- IT Deoxyribonucleotides
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(deoxyribodinucleotides, triphosphates, BODIPY
fluorophore-labeled; alternative dye-labeled primers,
ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for
automated DNA anal. and homogeneous amplification/detection assays)
- IT 187089-10-7DP, oligonucleotide primers labeled with
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(BODIPY 530/550; alternative dye-labeled primers,
ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for
automated DNA anal. and homogeneous amplification/detection assays)
- IT 150173-78-7DP, oligonucleotide primers labeled with
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(BODIPY 576/589; alternative dye-labeled primers,
ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for

automated DNA anal. and homogeneous amplification/detection assays)

IT 186961-29-5DP, oligonucleotide primers labeled with
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(BODIPY 589/616; alternative dye-labeled primers,
ribonucleotides, deoxyribonucleotides, and dideoxyribonucleotides for
automated DNA anal. and homogeneous amplification/detection assays)

IT 120718-52-7DP, TAMRA, oligonucleotide primers labeled with
138026-71-8DP, BODIPY, oligonucleotide primers labeled with
150152-69-5DP, BODIPY 581/591, oligonucleotide primers labeled
with 150173-72-1DP, oligonucleotide primers labeled with
150173-89-0DP, BODIPY 564/570, oligonucleotide primers labeled
with 165599-63-3DP, BODIPY 503/512, oligonucleotide primers
labeled with 174881-57-3DP, BODIPY 523/547, oligonucleotide
primers labeled with 186961-30-8DP, oligonucleotide primers labeled with
186961-31-9DP, oligonucleotide primers labeled with 186961-32-0DP
, oligonucleotide primers labeled with 186961-33-1DP,
oligonucleotide primers labeled with 186961-34-2DP, oligonucleotide
primers labeled with 186961-35-3DP, oligonucleotide primers labeled with
186961-36-4DP, oligonucleotide primers labeled with 186961-37-5DP,
oligonucleotide primers labeled with 186961-38-6DP, oligonucleotide
primers labeled with 186961-39-7DP, oligonucleotide primers labeled with
186961-40-0DP, oligonucleotide primers labeled with 186961-41-1DP,
oligonucleotide primers labeled with 186961-42-2DP, oligonucleotide
primers labeled with 186961-43-3DP, oligonucleotide primers labeled with
186961-44-4DP, oligonucleotide primers labeled with 186961-45-5DP,
oligonucleotide primers labeled with 186961-46-6DP, oligonucleotide
primers labeled with 186961-47-7DP, oligonucleotide primers labeled with
186961-48-8DP, oligonucleotide primers labeled with 186961-49-9DP,
oligonucleotide primers labeled with 186961-50-2DP, oligonucleotide
primers labeled with 186961-51-3DP, oligonucleotide primers labeled with
186961-52-4DP, oligonucleotide primers labeled with 186961-53-5DP,
oligonucleotide primers labeled with
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(alternative dye-labeled primers, ribonucleotides,
deoxyribonucleotides, and dideoxyribonucleotides for automated DNA
anal. and homogeneous amplification/detection assays)

L20 ANSWER 8 OF 8 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1990:548385 HCAPLUS

DOCUMENT NUMBER: 113:148385

TITLE: Incorporation of thionucleotides into nucleic acids
and oligonucleotides, and their application in nucleic
acid sequencing and hybridization assays using
fluorescence quenching

INVENTOR(S): Greulich, Karl Otto; Seidel, Claus; Wolfrum, Juergen;
Auer, Manfred; Gautel, Matthias; Goody, Roger; Labeit,
Siegfried

PATENT ASSIGNEE(S): Fed. Rep. Ger.

SOURCE: Ger. Offen., 7 pp.

CODEN: GWXXBX

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|------------------------|------|----------|-----------------|-------------|
| DE 3807975 | A1 | 19890928 | DE 1988-3807975 | 19880310 |
| DE 3807975 | C2 | 20020307 | | |
| US 6087101 | A | 20000711 | US 1997-990734 | 19971215 |
| PRIORITY APPLN. INFO.: | | | DE 1988-3807975 | A 19880310 |
| | | | US 1990-525038 | B1 19900518 |

AB A **fluorescent** dye which shows different quenching behaviors with the 4 nucleic acid bases is useful for detn. of base sequences in nucleic acids and oligonucleotides. Such **dyes** include **fluoresceins, rhodamines, coumarins, carbostyryls, and oxadiazoles**. The dye is attached via the S atom of a thionucleotide introduced e.g. by nick translation. Probes labeled in this manner are also useful in hybridization assays. Thus, a BglI-SalI fragment of the human immunodeficiency virus I pol reading frame was cloned in vector M13mp19, hybridized with GTAAAACGACGGCCA, and incubated (in 4 sep. reactions) with DNA polymerase Klenow fragment in the presence of TTP, dCTP, dGTP, dATP, and each of the 4 2',3'-deoxy-.alpha.-thionucleoside triphosphates S-labeled with 7-amino-N-(2-ethylaminocarbonyliodomethyl)-4-methylcoumarin as terminating nucleotides. The DNA segments were sepd. by PAGE in the presence of 8M urea in a "1-track-1-dye" method which used measurements of **fluorescence** lifetime of the dye to identify the base in each segment.

IT 19341-57-2


RL: PROC (Process)

(DNA incorporation of, for sequencing, **fluorescent** dye labeling in relation to)

RN 19341-57-2 HCAPLUS

CN Adenosine, 5'-(dihydrogen phosphorothioate) (8CI, 9CI) (CA INDEX NAME)

Absolute stereochemistry.



IC ICM C12Q001-68

ICS G01N033-50

ICA C07H021-04; C12P019-34; C07H021-00

CC 9-5 (Biochemical Methods)

ST **fluorescence** quenching nucleic acid sequencing; thionucleotide **fluorescent** label hybridization assay

IT Deoxyribonucleic acid sequences

(detn. of, thionucleotide incorporation and **fluorescence** quenching in relation to)

IT Nucleic acid hybridization

(thionucleotide incorporation and **fluorescent** dye labeling in relation to)

IT Nucleic acids

- RL: ANST (Analytical study)
(thionucleotide incorporation into, for sequencing, detection by
fluorescence quenching in relation to)
- IT Nucleotides, biological studies
RL: BIOL (Biological study)
(2',3'-dideoxy-, [.alpha.-thio]triphosphates, nucleic acid
incorporation of, for sequencing, detection by **fluorescence**
quenching in relation to)
- IT Dyes
(**fluorescent**, nucleic acid and oligonucleotide labeling with,
for hybridization assay, thionucleotide incorporation in relation to)
- IT Nucleotides, polymers
RL: ANST (Analytical study)
(oligo-, thionucleotides incorporation into, for sequencing and
hybridization assays, detection by **fluorescence** quenching in
relation to)
- IT Nucleotides, uses and miscellaneous
RL: USES (Uses)
(thio, nucleic acid incorporation of, **fluorescence** dye
labeling for hybridization assay in relation to)
- IT 15548-51-3 15867-02-4 **19341-57-2** 47151-76-8
RL: PROC (Process)
(DNA incorporation of, for sequencing, **fluorescent** dye
labeling in relation to)

10/049,164

=> d que

L45 802 SEA FILE=HCAPLUS ABB=ON PLU=ON (PROTEIN OR PROTEOM?) (3A) (ANAL Y? OR DETECT?) AND FLUORESC? AND (COMBINATOR? OR LIBRAR? OR PROBE)

L51 9 SEA FILE=HCAPLUS ABB=ON PLU=ON L45 AND MIXTUR? (5A) (PROTEIN OR PROTEOM?)

=> d ibib ab hitind 1-9

L51 ANSWER 1 OF 9 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2003:187663 HCAPLUS

TITLE: Peptidomics: A new approach to affinity protein microarrays

AUTHOR(S): Scrivener, Elaine; Barry, Richard; Platt, Albert; Calvert, Robert; Masih, George; Hextall, Patrick; Soloviev, Mikhail; Terrett, Jonathan

CORPORATE SOURCE: Oxford GlycoSciences, Abingdon, UK

SOURCE: Proteomics (2003), 3(2), 122-128

CODEN: PROTC7; ISSN: 1615-9853

PUBLISHER: Wiley-VCH Verlag GmbH & Co. KGaA

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Protein microarrays for diagnostic and **proteomic analyses** are being developed using a no. of different techniques for each of the steps required including immobilization methods, assay and detection systems. This is extremely different to the development of DNA microarrays which is now a well established technol. that has demonstrated the capabilities of transcriptomics to deliver validated differential transcripts. As mRNA and protein levels do not always correlate, protein microarrays would seem to be an obvious successor to DNA arrays. Unlike nucleic acids, however, protein targets are typically nonhomogeneous in physicochem. properties and affinity capture agents are often poorly characterised making the expts. difficult to perfect and reproduce. Moreover, running multiple affinity assays in parallel (multiplexing) is compromised by the heterogeneity of antibody affinities to their protein targets. In the peptidomic approach presented here the assayed **mixture of proteins** is enzymically digested prior to affinity capture to form a mixture of short peptides that are more similar in their physicochem. properties than intact proteins. These peptides can be predicted by in silico digestion of individual proteins, e.g. from protein databases allowing design of nonhomologous reagents for the screening of affinity agent **libraries**. The use of mass spectrometry (e.g. matrix-assisted laser desorption/ionization-time of flight mass spectrometry) for a direct confirmation of the identity of the species captured, provides a further advantage compared to the more usual method of detection in which **fluorescently** labeled captured species are scanned to give a spatially resolved image of the array.

CC 9 (Biochemical Methods)

L51 ANSWER 2 OF 9 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2003:113698 HCAPLUS

TITLE: Affinity **Analysis** of a **Protein** -Aptamer Complex Using Nonequilibrium Capillary Electrophoresis of Equilibrium Mixtures

AUTHOR(S): Berezovski, Maxim; Nutiu, Razvan; Li, Yingfu; Krylov,

CORPORATE SOURCE: Sergey N.
Department of Chemistry, York University, Toronto, ON,
M3J 1P3, Can.
SOURCE: Analytical Chemistry (2003), 75(6), 1382-1386
CODEN: ANCHAM; ISSN: 0003-2700
PUBLISHER: American Chemical Society
DOCUMENT TYPE: Journal
LANGUAGE: English

AB We propose a new method that allows the use of low-affinity aptamers as affinity **probes** in quant. **analyses** of **proteins**. The method is based on nonequil. capillary electrophoresis of the equil. **mixt.** (NECEEM) of a **protein** with its **fluorescently** labeled aptamer. In general, NECEEM of a protein with a **fluorescently** labeled aptamer generates an electropherogram with three characteristic features: two peaks and an exponential curve. Two peaks correspond to (i) the equil. amt. of free aptamer in the equil. **mixt.** and (ii) the amt. of the protein-aptamer complex that remains intact at the time of detection. The exponential part is ascribed to the complex decaying during sepn. under nonequil. conditions. Simple anal. of the three features in expts. with known concns. of the protein can be used for the detn. of the equil. dissocn. const., Kd, of the aptamer-protein complex. Similar anal. of the three features in the expt. with unknown concn. of the protein and known Kd value allows the detn. of the protein concn. In this proof-of-principle work, the NECEEM method was applied to the anal. of thrombin using a **fluorescein**-labeled aptamer under the conditions at which the protein-aptamer complex completely decayed during the sepn. We demonstrated that, despite the decay, as few as 4 .times. 10⁶ mols. of the **protein** could be **detected** with NECEEM without sacrificing the accuracy. This sensitivity is comparable with that reported by others for the aptamer-based equil. method. Thus, the proposed NECEEM-based method allows the use of aptamers for highly sensitive affinity **anal.** of **proteins** even when protein-aptamer complexes are unstable.

CC 9 (Biochemical Methods)

REFERENCE COUNT: 21 THERE ARE 21 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L51 ANSWER 3 OF 9 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2003:58701 HCAPLUS

DOCUMENT NUMBER: 138:119557

TITLE: Peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**

INVENTOR(S): Charych, Deborah; Beausoleil, Eric; Zuckermann, Ronald N.

PATENT ASSIGNEE(S): Chiron Corporation, USA

SOURCE: U.S. Pat. Appl. Publ., 32 pp., Cont.-in-part of U.S. Pat. Appl. 2002 55,125.
CODEN: USXXCO

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 3

PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|------------|------|------|-----------------|------|
|------------|------|------|-----------------|------|

| | | | | |
|-------|------|-------|-------|-------|
| ----- | ---- | ----- | ----- | ----- |
|-------|------|-------|-------|-------|

US 2003017508 A1 20030123 US 2002-190308 20020703
 US 2002055125 A1 20020509 US 2001-874091 20010604
 PRIORITY APPLN. INFO.: US 2000-209711P P 20000605
 US 2001-874091 A2 20010604

AB Provided are peptidomimetic protein-binding arrays, their manuf., use, and application. The protein-binding array elements of the invention include a peptidomimetic segment linked to a solid support via a stable anchor. The invention contemplates peptidomimetic array element **library** synthesis, distribution, and spotting of array elements onto solid planar substrates, labeling of complex **protein mixts.**, and the **anal.** of differential **protein** binding to the array. The invention also enables the enrichment or purifn., and subsequent sequencing or structural **anal.** of **proteins** that are identified as differential by the array screen. Kits including proteomic microarrays in accordance with the present invention are also provided. Slides were prepd. with a reflective aluminum coating that was further overcoated with a thin silicon dioxide dielec., followed by APTES. The Al/SiO2 substrate amplified the signal from Cy3/Cy5 tagged cDNA by approx. 10-40 fold relative to the corresponding glass substrate.

IC ICM G01N033-53

ICS G01N033-542; C12M001-34

NCL 435007900; 435287200

CC 9-1 (Biochemical Methods)

Section cross-reference(s): 3

ST peptidomimetic protein microarray mirrored substrate **proteomic analysis**; reflective aluminum silica APTES microarray

IT Proteins

RL: ARU (Analytical role, unclassified); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)

(A, biotinylated, immobilized on avidin coated slide; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)

IT Proteins

RL: ARU (Analytical role, unclassified); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)

(G, biotinylated, immobilized on avidin coated slide; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)

IT Organic compounds, uses

RL: DEV (Device component use); PRP (Properties); TEM (Technical or engineered material use); USES (Uses)

(aliph., C2-C100, as linker between anchoring and peptidomimetic segments on array; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)

IT Silanes

RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)

(amino, silicon dioxide modified with; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)

IT Thiols (organic), uses

RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)

(as anchoring agent; peptidomimetic protein-binding microarrays on

mirrored substrates for performing **proteomic analyses**

- IT Avidins
RL: ARU (Analytical role, unclassified); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)
(as coating on aluminum slides; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Polyoxyalkylenes, uses
RL: DEV (Device component use); PRP (Properties); TEM (Technical or engineered material use); USES (Uses)
(as linker or blocking agent on array; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Polymers, **analysis**
RL: ARU (Analytical role, unclassified); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)
(blocking agents; peptidomimetic **protein-binding** microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Electric insulators
(coatings, modified, on reflective metal on glass substrate; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT cDNA
RL: ANT (Analyte); ANST (Analytical study)
(**fluorescently**-labeled; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Proteins
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); BIOL (Biological study); USES (Uses)
(ligand-binding, peptidomimetic protein-binding, microarrays; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Peptides, uses
RL: DEV (Device component use); PRP (Properties); TEM (Technical or engineered material use); USES (Uses)
(oligopeptides, as linker between anchoring and peptidomimetic segments on array; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT DNA microarray technology
Glass substrates
Immobilization, molecular
Peptidomimetics
Protein microarray technology
(peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Proteome
RL: ANT (Analyte); BSU (Biological study, unclassified); PUR (Purification or recovery); ANST (Analytical study); BIOL (Biological study); PREP (Preparation)
(peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)

- IT mRNA
RL: ARG (Analytical reagent use); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)
(peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Proteins
RL: ANT (Analyte); BSU (Biological study, unclassified); PUR (Purification or recovery); ANST (Analytical study); BIOL (Biological study); PREP (Preparation)
(peptidomimetic; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Peptide library
(peptid; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Peptides, analysis
RL: ANT (Analyte); BSU (Biological study, unclassified); PUR (Purification or recovery); ANST (Analytical study); BIOL (Biological study); PREP (Preparation)
(peptoids, peptidomimetic; peptidomimetic **protein-binding** microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Metals, uses
RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)
(reflective, modified dielec. coating on; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT Amines, uses
RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)
(silyl, silicon dioxide modified with; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 6382-82-7
RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)
(aluminum slides coated with silicon oxide and layer of; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 541-59-3, Maleimide
RL: DEV (Device component use); PRP (Properties); RCT (Reactant); TEM (Technical or engineered material use); RACT (Reactant or reagent); USES (Uses)
(aminosilane functionalized with; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 89889-52-1
RL: DEV (Device component use); RCT (Reactant); TEM (Technical or engineered material use); RACT (Reactant or reagent); USES (Uses)
(aminosilane-coated aluminum slides coating with; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 7631-86-9, Silicon dioxide, uses
RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)
(aminosilane-modified, as coating on reflective metal on glass

- substrate; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 9013-20-1, Streptavidin 157885-16-0, Neutravidin
RL: ARU (Analytical role, unclassified); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)
(as coating on aluminum slides; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 919-30-2, APTES
RL: DEV (Device component use); RCT (Reactant); TEM (Technical or engineered material use); RACT (Reactant or reagent); USES (Uses)
(as coating on reflective aluminum slides; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 25322-68-3, Polyethylene oxide
RL: DEV (Device component use); PRP (Properties); TEM (Technical or engineered material use); USES (Uses)
(as linker or blocking agent on array; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 439084-54-5
RL: RCT (Reactant); RACT (Reactant or reagent)
(gold-coated microscope slides modification with; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 64987-85-5, SMCC
RL: RCT (Reactant); RACT (Reactant or reagent)
(in modification of gold-coated slides; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 58-85-5D, Biotin, conjugates with peptides or peptoids, immobilized on coated aluminum slides
RL: ARG (Analytical reagent use); DEV (Device component use); TEM (Technical or engineered material use); ANST (Analytical study); USES (Uses)
(peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 146368-14-1, Cy5 146368-16-3, Cy3
RL: RCT (Reactant); RACT (Reactant or reagent)
(proteins reaction with; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 7429-90-5, Aluminum, uses 7440-06-4, Platinum, uses 7440-16-6, Rhodium, uses 7440-32-6, Titanium, uses 7440-50-8, Copper, uses 7440-57-5, Gold, uses
RL: DEV (Device component use); TEM (Technical or engineered material use); USES (Uses)
(reflective, modified dielec. coating on; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)
- IT 221216-83-7
RL: PRP (Properties)
(unclaimed sequence; peptidomimetic protein-binding microarrays on mirrored substrates for performing **proteomic analyses**)

L51 ANSWER 4 OF 9 HCAPLUS COPYRIGHT 2003 ACS
 ACCESSION NUMBER: 2002:814272 HCAPLUS
 DOCUMENT NUMBER: 137:291270
 TITLE: Methods of **analysis** and labeling of
protein-protein interactions
 INVENTOR(S): Nollau, Peter; Mayer, Bruce J.
 PATENT ASSIGNEE(S): Children's Medical Center Corporation, USA
 SOURCE: PCT Int. Appl., 33 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|---------------|------|----------|-----------------|----------|
| WO 2002083846 | A2 | 20021024 | WO 2002-US11272 | 20020410 |
| WO 2002083846 | A3 | 20021212 | | |

W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN,
 CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,
 GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR,
 LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH,
 PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TN, TR, TT, TZ,
 UA, UG, US, UZ, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU,
 TJ, TM

RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH,
 CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR,
 BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG

PRIORITY APPLN. INFO.: US 2001-282748P P 20010410

AB We have discovered a new method to analyze and characterize complex cell signaling networks. The method is based on specific binding of protein-protein interaction modules to a single type of **protein** or a **mixture** of **proteins**. The method utilizes a no. of different protein-protein interaction domains as **probes** or sensors for the signaling state of the system under investigation. Glutathione-horseradish peroxidase conjugate was used to label fusion proteins of glutathione-S-transferase and protein-protein interaction domains. The labeled domains were applied to membranes on which lysates of 3T3 cells and v-abl transformed 3T3 cells were transferred. Signals were detected by chemiluminescence. Different patterns of tyrosine phosphorylated binding sites were detectable in different types of human leukemia when abl-SH2 or crk-SH2 were used as **probes** for the **detection** of **protein-protein** interactions.

IC ICM C12N

CC 9-2 (Biochemical Methods)

Section cross-reference(s): 6, 7, 14

IT Animal cell line

(3T3; methods of **anal.** and labeling of **protein**
 -protein interactions)

IT Protein motifs

(EH domain; methods of **anal.** and labeling of **protein**
 -protein interactions)

IT Protein motifs

(EVH1 domain; methods of **anal.** and labeling of
protein-protein interactions)

IT Proteins

RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);

ANST (Analytical study); BIOL (Biological study); USES (Uses)
(Grb-2, SH2 of, patterns of tyrosine phosphorylation in relation to;
methods of **anal.** and labeling of **protein-protein**
interactions)

IT Protein motifs
(LIM domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Protein motifs
(PDZ domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Protein motifs
(PTB domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Protein motifs
(RING finger; methods of **anal.** and labeling of
protein-protein interactions)

IT Protein motifs
(SAM domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Protein motifs
(SH2 domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT GTPase-activating protein
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
(SH2 of, patterns of tyrosine phosphorylation in relation to; methods
of **anal.** and labeling of **protein-protein**
interactions)

IT Protein motifs
(SH3 domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Protein motifs
(TPR domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Protein motifs
(WD40 repeat; methods of **anal.** and labeling of
protein-protein interactions)

IT Protein motifs
(WW domain; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Proteins
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
(abl, SH2 of, different patterns of tyrosine phosphorylated binding
sites in human leukemia in relation to; methods of **anal.** and
labeling of **protein-protein** interactions)

IT Biological materials
(**anal.** of; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Nucleic acids
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(analog, as label; methods of **anal.** and labeling of
protein-protein interactions)

IT Protein motifs
(ankyrin repeat; methods of **anal.** and labeling of
protein-protein interactions)

IT Protein motifs

(armadillo repeat; methods of **anal.** and labeling of **protein-protein** interactions)

IT Oligodeoxyribonucleotides
Oligonucleotides
Peptide nucleic acids
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(as label; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Membranes, nonbiological
(as support for **protein** immobilization; methods of
anal. and labeling of **protein-protein** interactions)

IT Plastics, reactions
RL: DEV (Device component use); RCT (Reactant); TEM (Technical or
engineered material use); RACT (Reactant or reagent); USES (Uses)
(as support for **protein** immobilization; methods of
anal. and labeling of **protein-protein** interactions)

IT Particles
(beads, as support for **protein** immobilization; methods of
anal. and labeling of **protein-protein** interactions)

IT Analysis
(biochem., multiplex binding assay; methods of **anal.** and
labeling of **protein-protein** interactions)

IT Proteins
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(conjugates, with oligonucleotides; methods of **anal.** and
labeling of **protein-protein** interactions)

IT Proteins
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
(crk, SH2 of, different patterns of tyrosine phosphorylated binding
sites in human leukemia in relation to; methods of **anal.** and
labeling of **protein-protein** interactions)

IT Oligonucleotides
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(derivs., thioesters, as label; methods of **anal.** and labeling
of **protein-protein** interactions)

IT **Fluorescent** substances
(fluor; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Immunoassay
(immunoblotting; methods of **anal.** and labeling of
protein-protein interactions)

IT **Proteins**
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
(labeled; methods of **anal.** and labeling of **protein**
-protein interactions)

IT Cell
Fluorescent substances
Human
Immobilization, molecular
Leukemia
Molecular association
PCR (polymerase chain reaction)
Protein motifs
Signal transduction, biological
(methods of **anal.** and labeling of **protein-protein**

- interactions)
- IT Antibodies
Fusion proteins (chimeric **proteins**)
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(methods of **anal.** and labeling of **protein-protein**
interactions)
- IT **Proteins**
RL: BSU (Biological study, unclassified); RCT (Reactant); BIOL (Biological
study); RACT (Reactant or reagent)
(methods of **anal.** and labeling of **protein-protein**
interactions)
- IT Phosphorylation, biological
(**protein**; methods of **anal.** and labeling of
protein-protein interactions)
- IT Gene, microbial
RL: BPN (Biosynthetic preparation); BSU (Biological study, unclassified);
BIOL (Biological study); PREP (Preparation)
(v-abl, 3T3 cells transformed with; methods of **anal.** and
labeling of **protein-protein** interactions)
- IT Protein motifs
(zinc finger; methods of **anal.** and labeling of
protein-protein interactions)
- IT 58-85-5, Biotin 50812-37-8, Glutathione-S-transferase
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(as label; methods of **anal.** and labeling of **protein**
-protein interactions)
- IT 9012-36-6D, Sepharose, conjugates with glutathione
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
(beads; methods of **anal.** and labeling of **protein**
-protein interactions)
- IT 50812-37-8D, Glutathione-S-transferase, fusion **proteins**
RL: ANT (Analyte); ARG (Analytical reagent use); BSU (Biological study,
unclassified); ANST (Analytical study); BIOL (Biological study); USES
(Uses)
(methods of **anal.** and labeling of **protein-protein**
interactions)
- IT 9003-99-0D, Peroxidase, conjugates with glutathione
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
ANST (Analytical study); BIOL (Biological study); USES (Uses)
(methods of **anal.** and labeling of **protein-protein**
interactions)
- IT 70-18-8DP, Glutathione, conjugates
RL: ARG (Analytical reagent use); BSU (Biological study, unclassified);
SPN (Synthetic preparation); ANST (Analytical study); BIOL (Biological
study); PREP (Preparation); USES (Uses)
(methods of **anal.** and labeling of **protein-protein**
interactions)
- IT 58-54-8, Ethacrynic acid 70-18-8, Glutathione, reactions
RL: RCT (Reactant); RACT (Reactant or reagent)
(methods of **anal.** and labeling of **protein-protein**
interactions)
- IT 60-18-4, Tyrosine, biological studies
RL: BSU (Biological study, unclassified); BIOL (Biological study)
(phosphorylated binding sites, different patterns of; methods of
anal. and labeling of **protein-protein** interactions)

L51 ANSWER 5 OF 9 HCAPLUS COPYRIGHT 2003 ACS
 ACCESSION NUMBER: 2002:696250 HCAPLUS
 DOCUMENT NUMBER: 137:228958
 TITLE: Protein profiling platform
 INVENTOR(S): Petricelli, Matthew
 PATENT ASSIGNEE(S): Activx Biosciences, Inc., USA
 SOURCE: PCT Int. Appl., 58 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|--|------|----------|-----------------|----------|
| WO 2002071066 | A1 | 20020912 | WO 2002-US6234 | 20020301 |
| W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG | | | | |
| US 2002182651 | A1 | 20021205 | US 2002-87602 | 20020301 |
| PRIORITY APPLN. INFO.: US 2001-273007P P 20010302 AB The invention concerns methods and compns. are described for analyzing complex protein mixts., such as proteomes, using activity-based probes. In particular, probes that specifically react with and bind to the active form of one or more target proteins are employed. Labeled peptides obtained from the labeled active target proteins can be used in screening and identification procedures, and can be related to the identity, presence, amt., or activity of active members of the desired target protein class. The methods and compns. described herein can be used, for example, to provide diagnostic information concerning pathogenic states, in identifying proteins that may act as therapeutic targets, and in drug discovery. IC ICM G01N033-53 ICS G01N033-543 CC 9-16 (Biochemical Methods) Section cross-reference(s): 1 IT Capillary electrophoresis Denaturants Diagnosis Diffusion Drug screening Electrophoresis Electrospray ionization mass spectrometry Fluorescent substances Gel electrophoresis HPLC Labels Liquid chromatography Mass spectrometry | | | | |

Protein degradation
 Test kits
 Washing
 (protein profiling platform)

REFERENCE COUNT: 3 THERE ARE 3 CITED REFERENCES AVAILABLE FOR THIS
 RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L51 ANSWER 6 OF 9 HCAPLUS COPYRIGHT 2003 ACS
 ACCESSION NUMBER: 2002:615942 HCAPLUS
 DOCUMENT NUMBER: 137:165832
 TITLE: Activity based **probe** analysis
 INVENTOR(S): Patricelli, Matthew P.
 PATENT ASSIGNEE(S): Activx Biosciences, Inc., USA
 SOURCE: PCT Int. Appl., 62 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|---------------|------|----------|-----------------|----------|
| WO 2002063271 | A2 | 20020815 | WO 2002-US3808 | 20020205 |
| WO 2002063271 | C1 | 20021024 | | |

W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN,
 CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,
 GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR,
 LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, OM, PH,
 PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TN, TR, TT, TZ,
 UA, UG, US, UZ, VN, YU, ZA, ZM, ZW, AM, AZ, BY, KG, KZ, MD, RU,
 TJ, TM
 RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH,
 CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR,
 BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG

PRIORITY APPLN. INFO.: US 2001-266687P P 20010205

OTHER SOURCE(S): MARPAT 137:165832

AB The invention concerns methods and compns. are described for
analyzing complex protein mixts. using
fluorescent activity-based probes. In particular,
probes that specifically react with and bind to the active form of
 one or more target proteins are employed. **Fluorescent** signals
 obtained from the labeled active target proteins can be related to the
 presence or amt. of active members of the desired target protein class.
 The methods and compns. described herein can be used, for example, to
 provide diagnostic information concerning pathogenic states, in
 identifying proteins that may act as therapeutic targets, and in drug
 discovery.

IC ICM G01N

CC 9-14 (Biochemical Methods)

Section cross-reference(s): 1, 14

ST protein sepn electrophoresis synthesis **fluorescent probe**
 drug screening

IT Electrophoresis

(SDS-PAGE; activity based **probe** anal.)

IT Capillary electrophoresis

Cyanine dyes

Diagnosis

Diffusion
Drug screening
Dyes
Electrophoresis apparatus
 Fluorescent substances
Fluorometry
Functional groups
Gel electrophoresis
Labels
Mass spectrometry
Pathogen
Separation
 (activity based **probe** anal.)
IT Receptors
 RL: ARG (Analytical reagent use); PRP (Properties); ANST (Analytical study); USES (Uses)
 (activity based **probe** anal.)
IT **Proteome**
 RL: BSU (Biological study, unclassified); PRP (Properties); BIOL (Biological study)
 (activity based **probe** anal.)
IT Functional groups
 (acylating; activity based **probe** anal.)
IT Functional groups
 (aldehyde; activity based **probe** anal.)
IT Functional groups
 (alkylating; activity based **probe** anal.)
IT Rare earth metals, uses
 RL: ARG (Analytical reagent use); PRP (Properties); ANST (Analytical study); USES (Uses)
 (cryptate derivs.; activity based **probe** anal.)
IT Functional groups
 (epoxide; activity based **probe** anal.)
IT Functional groups
 (ketone; activity based **probe** anal.)
IT Dyes
 (metal chelate; activity based **probe** anal.)
IT **Proteins**
 RL: ANT (Analyte); BSU (Biological study, unclassified); PRP (Properties); ANST (Analytical study); BIOL (Biological study)
 (**mixt.**; activity based **probe** anal.)
IT Dyes
 (naphthylamine; activity based **probe** anal.)
IT Reagents
 RL: NUU (Other use, unclassified); USES (Uses)
 (noncovalent; activity based **probe** anal.)
IT Functional groups
 (phosphoryl; activity based **probe** anal.)
IT Functional groups
 (sulfonyl; activity based **probe** anal.)
IT 13558-31-1 98181-63-6
 RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
 (activity based **probe** anal.)
IT 91-64-5, Coumarin 92-83-1, Xanthene 7440-18-8D, Ruthenium, chelates 7440-27-9D, Terbium, chelates 7440-52-0D, Erbium, chelates 25168-10-9, Naphthylamine 138026-71-8, BODIPY
 RL: ARG (Analytical reagent use); PRP (Properties); ANST (Analytical

study); USES (Uses)
(activity based **probe** anal.)

IT 446828-34-8P 446828-36-0P 446850-41-5P 446850-43-7P 446850-45-9P
446850-47-1P
RL: ARG (Analytical reagent use); RCT (Reactant); SPN (Synthetic
preparation); ANST (Analytical study); PREP (Preparation); RACT (Reactant
or reagent); USES (Uses)
(activity based **probe** anal.)

IT 189200-71-3DP, Rhodamine green, reaction with adenosine derivs.
446833-62-1P 446833-64-3P 446850-50-6P 446850-53-9P 446850-55-1DP,
reaction with rhodamine green 446850-58-4P 446850-61-9P 446850-64-2P
446850-67-5P 446850-69-7DP, reaction with rhodamine green 446850-71-1P
446850-73-3P 446850-76-6P 446850-79-9P 446850-81-3P
RL: ARG (Analytical reagent use); SPN (Synthetic preparation); ANST
(Analytical study); PREP (Preparation); USES (Uses)
(activity based **probe** anal.)

IT 112-47-0, 1,10-Decanediol 112-60-7
RL: RCT (Reactant); RACT (Reactant or reagent)
(activity based **probe** anal.)

IT 134179-40-1P 338964-01-5P 338964-02-6P 338964-03-7P 338964-04-8P
338964-05-9P 338964-06-0P
RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT
(Reactant or reagent)
(activity based **probe** anal.)

L51 ANSWER 7 OF 9 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2002:564598 HCAPLUS

TITLE: Proteomic profiling of mechanistically distinct enzyme
classes using a common chemotype

AUTHOR(S): Adam, Gregory C.; Sorensen, Erik J.; Cravatt, Benjamin
F.

CORPORATE SOURCE: The Skaggs Institute for Chemical Biology and
Department of Chemistry, The Scripps Research
Institute, La Jolla, CA, 92037, USA

SOURCE: Nature Biotechnology (2002), 20(8), 805-809
CODEN: NABIF9; ISSN: 1087-0156

PUBLISHER: Nature Publishing Group

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Proteomics research requires methods to characterize the expression and
function of **proteins** in complex **mixts**. Toward this
end, chem. **probes** that incorporate known affinity labeling
agents have facilitated the activity-based profiling of certain enzyme
families. To accelerate the discovery of proteomics **probes** for
enzyme classes lacking cognate affinity labels, we describe here a
combinatorial strategy. Members of a **probe**
library bearing a sulfonate ester chemotype were screened against
complex proteomes for activity-dependent protein reactivity, resulting in
the labeling of at least six mechanistically distinct enzyme classes.
Surprisingly, none of these enzymes represented targets of previously
described proteomics **probes**. The sulfonate **library**
was used to identify an omega-class glutathione S-transferase whose
activity was upregulated in invasive human breast cancer lines. These
results indicate that activity-based **probes** compatible with
whole-**proteome anal.** can be developed for numerous
enzyme classes and applied to identify enzymes assocd. with discrete
pathol. states.

CC 7-1 (Enzymes)
Section cross-reference(s): 6, 13, 9, 14
ST proteome protein enzyme profile **fluorescent probe**
IT **Fluorescent** indicators
Human
(proteomic profiling of mechanistically distinct enzyme classes using a
common chemotype)
REFERENCE COUNT: 29 THERE ARE 29 CITED REFERENCES AVAILABLE FOR THIS
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L51 ANSWER 8 OF 9 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1993:142804 HCAPLUS

DOCUMENT NUMBER: 118:142804

TITLE: **Analysis** of a recombinant **protein**
preparation on physical homogeneity and state of
aggregation

AUTHOR(S): Brochon, J. C.; Tauc, P.; Merola, F.; Schoot, B. M.

CORPORATE SOURCE: LURE, Cent. Univ. Paris-Sud, Orsay, F91405, Fr.

SOURCE: Analytical Chemistry (1993), 65(8), 1028-34

CODEN: ANCHAM; ISSN: 0003-2700

DOCUMENT TYPE: Journal

LANGUAGE: English

AB The homogeneity of a recombinant protein prepn. and the state of
aggregation were studied by time-resolved polarized fluorometry. From the
rotational correlation time .theta. distribution pattern, the state of
aggregation of a protein in soln. is deduced. This distribution is detd.
from a 2-dimensional (.tau., .theta.) fit using the max. entropy method of
data anal. and where .tau. values are the **fluorescence**
lifetimes. An anal. procedure is developed that is validated by
measurements on a **mixture** of 2 **proteins** having different
mol. wts. and contg. a single tryptophan residue per polypeptide chain:
recombinant human interferon .gamma., r-hu IFN .gamma. (RU42369) of Mr
17,000 for the monomer and recombinant mutant W201Y lac operon repressor
tetramer, Mr 152,400. By spiking a soln. of 1 mg/mL r-hu IFN .gamma. with
the W201Y lac operon repressor, the lower level of detection of higher
mol. wt. component is found to be 5% in **intrinsic fluorescent**
probe concn. In this study it was found that (1) purified r-hu
IFN .gamma. in soln. after lyophilization is a dimeric mol., without
indication of phys. heterogeneity and without high-mol.-wt. aggregates;
(2) heat treatment of lyophilized r-hu IFN .gamma., 14 days at 40.degree.,
results in the formation of a detectable amt. of higher-mol.-wt. material;
(3) the dissocn. of the dimer r-hu IFN .gamma. on diln. was not detectable
after diln. to 0.01 mg/mL (0.57 .mu.M). Taking advantage of the great
sensitivity of a **fluorescence** technique and of the capabilities
of the data anal. MEM, this new procedure can be widely used to detect a
high-mol.-wt. protein contaminant (aggregates) in a homogeneous protein
soln.

CC 9-5 (Biochemical Methods)

Section cross-reference(s): 16

ST recombinant protein aggregation homogeneity detn; fluorometry recombinant
protein analysis

IT Homogeneity

Molecular association

(of recombinant **proteins**, fluorometric anal. of)

L51 ANSWER 9 OF 9 HCAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1991:531404 HCAPLUS

DOCUMENT NUMBER: 115:131404
 TITLE: Luminescent enzyme assay methods and kits for detecting a substance using enzymically-induced decomposition of dioxetanes
 INVENTOR(S): Bronstein, Irena Y.
 PATENT ASSIGNEE(S): Tropix, Inc., USA
 SOURCE: U.S., 30 pp. Cont.-in-part of U.S. Ser. No. 265,406, abandoned.
 CODEN: USXXAM
 DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 16
 PATENT INFORMATION:

| PATENT NO. | KIND | DATE | APPLICATION NO. | DATE |
|--|------|----------|-----------------|-------------|
| US 4978614 | A | 19901218 | US 1989-382125 | 19890720 |
| CA 2033331 | AA | 19910121 | CA 1990-2033331 | 19900717 |
| WO 9101492 | A1 | 19910207 | WO 1990-US3920 | 19900717 |
| W: CA, JP | | | | |
| RW: AT, BE, CH, DE, DK, ES, FR, GB, IT, LU, NL, SE | | | | |
| EP 435998 | A1 | 19910710 | EP 1990-911242 | 19900717 |
| EP 435998 | B1 | 19991110 | | |
| R: DE, FR, GB, IT | | | | |
| EP 859062 | A2 | 19980819 | EP 1998-101157 | 19900717 |
| R: DE, FR, GB, IT | | | | |
| EP 907082 | A2 | 19990407 | EP 1998-122597 | 19900717 |
| R: DE, FR, GB, IT | | | | |
| US 5220005 | A | 19930615 | US 1990-574787 | 19900830 |
| JP 04124186 | A2 | 19920424 | JP 1990-239765 | 19900910 |
| JP 2999810 | B2 | 20000117 | | |
| US 5605795 | A | 19970225 | US 1994-255795 | 19940607 |
| US 5856522 | A | 19990105 | US 1997-882330 | 19970625 |
| US 36536 | E | 20000125 | US 1997-958342 | 19971027 |
| PRIORITY APPLN. INFO.: | | | | |
| | | | US 1988-265406 | B2 19881026 |
| | | | US 1983-889823 | A2 19830724 |
| | | | US 1986-889823 | A2 19860724 |
| | | | US 1989-382125 | A 19890720 |
| | | | EP 1990-911242 | A3 19900717 |
| | | | WO 1990-US3920 | W 19900717 |
| | | | US 1990-574787 | A3 19900830 |
| | | | US 1992-990800 | B1 19921215 |
| | | | US 1995-433996 | A1 19950504 |

OTHER SOURCE(S): MARPAT 115:131404

AB Luminescent enzyme assays and kits are described that use an enzyme and dioxetane I (T = cycloalkyl or polycycloalkyl group bonded by a spiro linkage; Y = **fluorescent** chromophore; X = H, alkyl, aryl, aralkyl, heteroalkyl, heteroaryl, cycloalkyl, cycloheteroalkyl, or enzyme-cleavable group; Z = H, enzyme-cleavable group; .gtoreq.1 of X and Z must be an enzyme-cleavable group) so that the enzyme cleaves the enzyme-cleavable group from the dioxetane to form a neg.-charged substituent bonded to the dioxetane. The neg.-charged substituent causes the dioxetane to decomp. to form a luminescent substance comprising group Y. Human chorionic gonadotropin (hCG) was detd. in blood and urine by chemiluminescence immunoassay using 3-(2'-spiroadamantane)-4-methoxy-4-(3''-phosphoryloxy)phenyl-1,2-dioxetane, di-Na salt (AMPPD) as a substrate, alk. phosphatase conjugated with anti-hCG antibody, and beads

coated with anti-hCG. Samples were read in a luminometer. The assay was >10 times as sensitive as a colorimetric assay using p-nitrophenylphosphoric acid as a substrate. In an assay for alk. phosphatase using AMPPD, the min. detectable concn. of alk. phosphatase was 1.67 .times. 10⁻¹⁵M.

IC ICM G01N021-76

ICS G01N033-53

NCL 435021000

CC 9-5 (Biochemical Methods)

Section cross-reference(s): 2, 7

IT Antibodies

Antigens

Proteins, analysis

RL: ANT (Analyte); ANST (Analytical study)

(detn. of, by luminescence enzyme immunoassay, dioxetane substrate in)

IT Enzymes

RL: ANST (Analytical study)

(conjugates, with antibody or hybridization **probe**-linkable

agent, in enzyme luminescence assays with dioxetane substrate)

IT Albumins, compounds

RL: ANST (Analytical study)

(conjugates, with **fluorescein**, alk. phosphatase luminescence

assay with dioxetane substrate response to)

IT 2321-07-5D, **mixts.** with polymeric quaternary ammonium salts

9017-80-5, Poly(vinylbenzyltrimethylammonium chloride) 135781-06-5

135781-07-6 135781-08-7

RL: ANST (Analytical study)

(**protein** luminescence enzyme assay with dioxetane substrate enhancement with)